





Comfort  
Zone  
Technology

### Takes Your Mind Off The Pain

Pressing Nova Safety Lancets' Comfort Zone of 8 raised dots against the skin sends a message of comfort to the brain. Then, when you click the release button, the lancet penetrates your skin, but the dominant Comfort Zone sensation means the pain 'message' to your brain is greatly reduced. Distracting the brain in this way makes sampling virtually pain-free.\*

\* Data on file.

## Nova Safety Lancets Offer Two Sizes to Meet Your Point-of-Care Needs

Lancet Type	Number In Case	Spec	Order Number
 Nova Safety Lancet Comfort	4000	28G Depth 1.8mm Avg. Bld. vol 5µL-30µL	50938
 Nova Safety Lancet Normal	4000	23G Depth 1.8mm Avg. Bld. vol 30µL-75µL	50939

For more information please call Nova Biomedical  
1-800-458-5813, or visit [www.novabiomedical.com](http://www.novabiomedical.com)

## nova Safety Lancets



Nova Safety Lancet  
Comfort  
28G



Nova Safety Lancet  
Normal  
23G

Nova Safety Lancets Offer  
Two Sizes to Meet Your  
Point-of-Care Needs

**nova**  
biomedical  
novabiomedical.com

200 Prospect Street, Waltham, MA 02454-9141 U.S.A.  
TEL: (781) 894-0800 (800) 458-5813 FAX: (781) 894-5915  
Int'l FAX: (781) 899-0417

Comfort Zone Technology is a registered trademark of Owen Mumford, Ltd.

No. 306 12-22-11

## How to Make Blood Sampling Less of a Pain

## nova Safety Lancets

Single-use  
Comfort Zone Technology®  
from Nova Biomedical®



**nova**  
biomedical

**nova**  
**Safety Lancets**

Easy, comfortable blood sampling – now available from Nova Biomedical®


With Nova Safety Lancets there's now a convenient new way to test your blood glucose levels that's a lot less painful. That's because Nova Safety Lancets' unique Comfort Zone Technology® actually distracts your brain from feeling the prick of the needle.

**Safe and Hygienic**

As a single-use safety lancet, you will reduce your risk of infection by using a fresh Nova Safety Lancet every time. There's no need to load the lancet, and you never see it because it remains hidden inside the case. So there's also little chance of accidental injury before, or after, the test.




Nova Safety Lancet  
Comfort  
28G




Nova Safety Lancet  
Normal  
23G


**Nova's Safety Lancets Offer Two Sizes to Meet Your Patient Needs**




**TWIST**



**CLICK**



**TEST**



**DISPOSE**

**Using Nova Safety Lancets couldn't be simpler**

Simply twist off the protective sterile cap, press the lancet against your finger, click the release button, test with any meter device and dispose. It's that easy.

