



Critical Care Animal Testing Made Simple

Introducing a new generation veterinary critical care analyzer that combines the revolutionary micro-electronics of the consumer world with Nova Biomedical's innovative MicroSensor Card™ for a simple, compact, fast, and economical analyzer.

10-Test Stat Menu

pH, PCO_2 , PO_2 , Na, K, Cl, iCa
Glu, Lac, Hct

100 Microliter Blood Sample

Results in 1 minute

Clot Block™ Protection

Unique probe and flow path protects the analyzer from blood clots.

Total Quality Monitoring

- Automated Quality Control verifies analyzer accuracy and precision on assayed liquid controls.
- Supplemental Quality Monitoring electronically verifies the performance of all analytical components on every analysis.



Maintenance-Free
MicroSensor Card
2.5 in x 2 in (6.35 cm x 5.0 cm)

Keep Profits In-House for Blood Testing



**Essential 10-test menu
optimizes blood testing
revenues**

**Open data interface and
optional data management
system optimize billing**

**Maintenance-free system
keeps clinical staff focused
on patient care**

Compact, Portable Size

Stat Profile Prime VET is one of the smallest and lightest critical care analyzers. It is so compact it can be located anywhere in the veterinary hospital or operated on a mobile cart.

